


# STD Assays

## STD Kit



Detection and differentiation of five different STD pathogens on one strip.

Each PCR-Mix can be used as a stand-alone assay or in any combination, as required.

An inbuilt control (Amplification control) ensures that DNA has been extracted and amplified from the starting test material.


Detection of 17 high and 6 low risk HPV genotypes

Detection of *Chlamydia trachomatis* with deleted or mutated variant in the cryptic plasmid

Specimens material: cervical or urethral swab, urine and samples coated with paraffin wax

CE 0123 certified

## Chlamydia trachomatis Kit



**Detection of HPV and HSV also available as single testkits!**

Order-No.:

STD	RDB2110(X)	12(60) tests
C.t.	RDB2117(X)	12(60) tests
HPV screening	RDB2082(X)	12(60) tests
HPV typing	RDB2080(X)	12(60) tests
HSV1/2	RDB2130(X)	12(60) tests

**Ready to use reagents**